

Expression of cellulases in plants...

Rationale: For an initial set of constructs, fuse E2 and E3 genes of Thermoplasma fusca to MAC promoter and NOS terminator in pCGN1578 constructs similar to Dennis-321 (MnF) or 263 (LiP) but lacking any signal sequence. Based on published information, neither protein is heavily glycosylated (E3 is functional without glycosylation), so cytosolic localization should be compatible with function.

9/18/95

streaked from agar slants to LBamp plates
 D644 gave only 3 colonies (DHSa + pSZ10)
 D651 " many " (" pSZ15)

patch 2 of each for mini-preps.

1/19/96 - prep 1/4 plate each prep
 using "Vince Schultz" protocol
 (NH₄NO₃ mtn., G10H mtn.)

- resuspended in 50 μ l, load 2 μ l

From maps...

E2 - 348-1578

Sma I 2 sites
 Sal I
 Nar I 2 sites
 Eco RI
 Sph I
 Sac II
 Bst E II 2 sites
 Bp M I
 Pvu II

E3 - 689-2363

Kpn I
 Bp M I "
 Bst E II "
 Sal I "
 Nar I "
 - Sac II
 Nco I
 Bst X I
 Sac I "
 Hpa I
 Sma I

EXHIBIT

B

Thomas Ziegelhoffer

E3 construct - PCR amplification

Overall strategy:

- 1) generate XbaI @ 5' end, EcoRI @ 3' end (PCR)
- 2) subclone into pUC19 as XbaI-EcoRI fragment
- 3) replace ~1.4 kb KpnI-SmaI of resulting pUC-E3 clone with KpnI-SmaI of pJZ10 (from strain D(44) → sequence using universal primers)
- 4) clone XbaI-EcoRI fragment of fully wt. pUC-E3 into 253-1 vector (3 way ligation - XbaI-EcoRI-HindIII)
- 5) clone expression cassette (MAC-E3-MASTER) into pCGN1578 as BamHI-HindIII fragment

Custom primers from UWBC (see data sheets)

XbaE3 (DSN# 5500) - designed to create an XbaI site and ATG initiation codon at the normal leader cleavage site, generating the sequence NH₂-MET-ALA-GLY-Cys-Ser-...
 25mer of ~16225 MW, resusp. to give 210 μ l solution (348 μ g in 100 μ l)

RIE3 (DSN# 5499) - designed to create an EcoRI site immediately after the stop codon of E3.
 23mer of ~14927 MW, resusp. to give 205 μ l solution (307 μ g in 100 μ l)

Set up 1st 6 of buffer optimization set


10 μ l	10x buffer (#s 1-6)	}	combine
10 μ l	4 dNTPs (2.5mM each)		
1 μ l	template (1/1000 dilution of pJZ10 prep)		
1 μ l	10 μ M XbaE3		
1 μ l	10 μ M RIE3		
7.5 μ l	dH ₂ O		
0.25 μ l	Taq pol	470 μ l dH ₂ O	

Method 1 change prog 31 to 48°C annealing temp.


Theresa Fisher

E2 cloning - strategy

E2 gene (in pUC19)




253-1 expression cassette (in pUC118, I think...)



Advantage: RI site available immediately downstream of 3' end of E2.
 Problem: Need to avoid using EcoRI (2 sites in 253-1, 1 site in E2 gene)

final step in expr. cass. construction



- vector backbone is 253-1, cut with XbaI & HindIII
- Xba-Sph fragment (~375 bp) is PCR-derived and contains the N-terminal end of E2 with new N-terminal Met (Met-Asn-Asp-Ser... etc). This fragment will also be subcloned into pUC19 for sequencing to verify absence of errors
- Sph-Hind fragment (1.6 kb) is derived from intermediate clone (~1.0 kb RI fragment of E2 cloned into RI site of 253-1 [~2.5 kb RI fragment excised])

Thomas J. ...

E2 & E3 - more pCGN cloning

resolve mini-preps of colonies on 0.5% gel

was plasmid! (with one exception!)

#19 is coincidentally the only isolate which grew significantly upon patching to G20 (others are reversion background)

cut 1 μ l of 19 DNA
 0.5 μ l cox multi
 0.25 μ l SmaI
 1 μ l #19
 7.25 μ l dH_2O

Result:

LiP - 2.0, 2.4, 2.9, 3.2, 7.0
 total = 17.5 kb


#19 - 1.4, 2.9, 3.0, 3.2, 7.0
 total = 17.5 kb

Q: is pCGN1578 map correct?

do more digests but first,


transform comp LBA4404 w/ 5 μ l
 express 24 hr @ 28°C
 plate on YEP 20G @ 28°C

Thomas J. Leggett



migrator distance	number	elimer
4.3 kb	3.6	3.1
7.0 kb	3.2	2.7
17.5 kb	2.4	1.7
2.7 kb	—	3.3

see reverse for graph --



Move E2 & E3 in pCGN1578 ---

digest 318-S with SmaI

digest 260-6, 318-S, #19 with EcoRI (multi)

all should have 1.0, 4.6, 7.3

260, 318 should also have ~2.4 kb fragment, ~2.0 kb
E2 " " " ~0.95, 1.45, ~2.0 kb


Result: PERFECT!

318-S - fragments of
~2.4, 2.7, 3.2, 9.5 (180000)
(with SmaI)

others are as predicted

CONCLUSION:

- map of pCGN1578
is missing a SmaI
site @ ~1600-1750
- E2 construct (#19)
is correct!



THEREFORE,

- bona fide gent^R transformants are recovered at very low frequency
- reversion frequency (spontaneous DH5a gent^R) is significant ($\sim 10^{-6}$ - 10^{-7}) but these are much weaker growers than gent^R transformants

X Thomas J. [Signature]

X

E2 & E3 - pCGN constructs

using ligation mixtures described on p. 20, retransform comp DH5 α and plate everything...

- E2 1 ml on 3 plates
- E3 2 ml on 6 plates

Result: E2 15 lg. colonies
E3 5 " " (3 lg., 2 weak, actually)

Note: control transformation with 260-6 shows that lg. colonies correspond to true transformants (efficiency of transf. with 260-6 was comparable to pMD4.21, a pBR322 derivative)

start 1.5 ml LB 20G cultures for mini-prep (3 E2 clones, 6 E3 clones)

harvest after ~8 hr growth @ 33°C
E2-1-3, E3-1-3 are good E3-5,6 are weak
and E3-4 is dead

3/24 - mini-prep on cell pellets

digest with EcoRI in "H" buffer

all E2 isolates look good
all E3 " " bad

- poor growers (on 20G) #5,6 appear to have no plasmid (as expected from previous results)
- #1,2 are indistinguishable from original 260-6 clone and probably represent background of single cut vector

Thomas J. [signature]

pCGN-E3, the continuing saga...

Since no "real" E3 transformants were recovered,
try again

- use vector purified ^{as} previously (p20)
- use insert " " " (p15)

Ligation:	old	new
vector	4	6.5
insert	2	2
10x buffer	1	1
T4 ligase	0.5	0.5
dH ₂ O	2.5	—

transformants:

test by EcoRI digestion of miniprep DNA


#s 2, 7, 10, 16 are correct!
with the exception of #1
(?!), all others are
260-6 recovered

transform #2 into camp type
(2 μ l DNA + 100 μ l cells)

approx 500 transformants
obtained (YEP 20G)

start culture for plant
transformation

Thomas Ziegler



E2 - Western blot (std. dilutions, callus, leaves)

conc. of E2 is 1mg/ml

- dilute 1:100 in 0.1mg/ml BSA
- load 1, 2, 5, 10 μ l of diluted enzyme (10, 20, 50, 100 ng)

grind tissue (Tobacco transformation) in 1x sample

use 1:50 dilution of anti E2 Ab (Diana Irwin's recommendation)

very good signal @ 10ng E2

- approximately 20 μ l of tissue (20mg)
- assuming ~1% as soluble protein = 200ng total
- load ~ $\frac{1}{4}$ th of total (=50ng)
- signal looks like ~ 3-10ng (ie. - a high % of total!)

VERY PROMISING! (but is it active?...)

* plant-expressed E2 is ~ same MW as purified enzyme

E2 - western blot

as on p27, dilute E2 stock in dH_2O / BSA ^{0.1 mg/ml}
 1, 2, 5, 10 μl of 1:100 dilution (10-100 ng)

test selected extracts from p29

Result: 6/3

lane	
1	MW marker
2	10 μl E2 d.i.
3	5 " "
4	2 " "
5	1 " "
6	CT 1
7	" 2
8	" 8
9	" 10
10	" 13
11	" 25
12	" 30
13	" 38
14	" 39
15	" 30 (1:10 dilution in dH_2O + BSA) 0.2 μl

10% gel

1 μl
each

none detect.
 none detect.
 < 1 ng
 << 1 ng
 none detect.
 ~ 2 ng
 ~ 3 ng
 ~ 2 ng
 none detect.
 < 1 ng

Result: Western blot agrees (for the most part) with activity measurements, with some notable exceptions:

- ① based on CMC, CT13 should have yielded protein
- ② based on MUCB, CMC CT39 " " " "
- ③ CT25, 38 show good protein yield, despite low MUCB activity

- based on MUCB assay, expected ~ 500 ng E2
- observed less than 10 ng, therefore activity is a very poor measure of expression
- only about 0.1% total extracted protein is E2

Thomas Ziegler

Update on plant analysis - W38/MAC-E2 (TZA1)

Analysis to date: ① reducing sugar assay using CMC as substrate, ② MUCB fluorometric assay substrate ③ Western blot using anti-E2 antiserum.

Results: ① based on Western blotting, E2 expression is @ ~ 0.1% total extracted protein for CT25, CT30, CT38 and is barely detectable in CT8, CT10. Other plants tested "contained undetectable amounts of E2". ② Although there is some correlation between CMC or MUCB activity and the presence of E2 (as determined by Western), the values obtained are orders of magnitude higher than would be expected, suggesting synergistic activity with endogenous cellular activities. ③ Using pure E2 provided by D. Irwin, approximate K_m , V_{max} for MUCB are $\sim 9 \mu M$ and ~ 0.12 mole/min mole, respectively.

Conclusions: A better assay method is needed.

① CMC/reducing sugar assay is probably not practical for lg. # of samples ② MUCB assay could be improved by (A) extracting juice directly (no buffer) followed by (B) heating @ $60^\circ C$ for 30 minutes. Lower conc. of substrate (C) in the $20-100 \mu M$ range and the addition of $1 mM$ glucosyl-1,5-lactone to inhibit β -glucosidase activities may also help. Could try filter adsorption of heat-treated sample to further clean up (Millipore Ultrafree-DEAE_{mc}). ③ For the time being, screen plants by Western blot until assay is improved.

Thomas Feghly

E2 Western blot ^{TZA1} (W38/MAC_p-E2)

grind ~ 20 mg leaf tissue in 40 μ l ~~buffer~~ ^{0.1% SDS also} buffer
(50 mM NaOAc, pH 5.5, 100 mM KCl). Load 2 μ l
of extract / lane. save remainder for juice
extraction / MUCB assay (freeze leaf t. @ -80°C)

1	(control)	1	CT 85	★
2	1 μ l of E2 dilution (10 μ g)	2	CT 86	★
3	CT 6	3	CT 87	★
4	CT 9	4	1 μ l of E2 dilution	
5	CT 11	5	CT 88	★
✓ 6	CT 15	6	CT 89	
7	CT 16	7	CT 90	★
✓ 8	CT 17	✓ 8	CT 91	
✓ 9	CT 18	9	CT 92	
10	CT 20 ★	10	CT 94	★
11	CT 22	11	CT 95	
✓ 12	CT 26	12	CT 96	★
13	CT 36 ★	✓ 13	CT 97	
14	CT 40 ★	14	CT 99	★
15	CT 84	15	CT 100	★

10% gel → transfer → Western (anti-E2 solution, used 2x prev.)

RESULT: Those marked with a check show little or
no E2 protein. ★ indicates doublet (modification?)
Extracts are ~ 1 mg/ml. Estimate that CT 96, CT 99
are ~ 2 μ g phytoe in 2 μ g total protein.

Thomas J. J. J.

~~XXXXXXXXXX~~

E3 western blot (tobacco transformants)

Grind ~40mg leaf tissue in 80 μ l "buffer D" (see p. 44), spin down pellets, recover sup (done by BK on 7/8/96).

prepare samples: 20 μ l extract + 5 μ l 5x buffer

Lane#	gel 1	E3? ↓	gel 2	E3? ↓	E3 diluted into dH ₂ O + 0.1mg/ml BSA
1	MW marker				
2	10 μ l E3 1:100		8 μ l CT60 sample	-	
3	5 μ l " "		61	++	10% gels
4	2 μ l " "		63	+/-	
5	1 μ l " "		64	-	
6	8 μ l CT47 sample		66	-	1:200 of E3 anti-serum
7	" 48	?	68	+	
8	" 49	+	69	-	
9	" 51	+	71	+/-	
10	" 52	+/-	72	+	
11	" 53	+/-	75	+++	
12	" 54	++	76	+/-	
13	" 55	++	77	+/-	
14	" 57	++	78	"lost" (false lane)	
15	" 59	-	79	++	

Result: all are significantly less than 10 μ g control. CT75 appears to be highest (maybe ~5 μ g). A cross-reacting band @ ~30 kDa was present in all lanes.

Thomas J. Dwyer

4/18

10 20 50 100

$$\underline{21C} \quad 2714C$$

6050

$\frac{d}{dt} \left(\frac{1}{\rho} \right) = - \frac{1}{\rho^2} \frac{d\rho}{dt}$

E-2-2

← E2

6/17

CT-1 control

6 9 11 15 16 12 18 20 22 26 36 40 84 85 86 87 88 89 90 91 92 94 95 96 97 99 100

E2 - 1

 $\leftarrow E2$

ng E3

100 50 20 0 47 48 49 51 52 53 54 55 57 59

60 61 63 64 66 68 69 71 72 75 76 77 78 79

$$E \rightarrow$$

4 E3